

Quantist™ Release Notes

P/N 805-0004

Revision B, October 2023

Copyright © 2023 Bio-Techne. All rights reserved.

Version 1.0.1

Change History

Change Type	Notes
-------------	-------

Improvement	Patch Release
Bug fix	

Fixed Issues:

- Loading Luminex CSV File with standard dilutions all containing the same concentration will cause application to become unresponsive (BTQ-572)
- Support installation with non English Windows 10 OS (BTQ-576)

Version 1.0.0

Change History

Change Type	Notes
-------------	-------

Feature	Initial Release
---------	-----------------

Known Issues

- Well Assignments View can crash application if assignments are performed quickly back-to-back (BTQ-569)
- Loading Luminex CSV File with standard dilutions all containing the same concentration will cause application to become unresponsive (BTQ-572)

